EPflex

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ENGLISH

Note Carefully read and understand the entire Instruction for use prior to using the guidewire. Failure to do so could result in complications.

Intended purpose

The Guidewires are intended to facilitate the introduction and placement of catheters or other interventional or diagnostic devices within the central circulatory system, the heart chambers and the coronary vasculature during diagnostic or interventional procedures.

Indications

The guidewires are indicated for use during the following clinical procedures:

- Percutaneous Transluminal Coronary Angioplasty (PTCA)
- Transcatheter Aortic Valve Implantation (TAVI)
- Placement of septal occlude
- Endovascular aneurysm Repair (EVAR)
- Placement of vascular catheters e.g. CVC, PICC, CVWH

ATTENTION: Only guidewires labelled 'TAVI application' may be used for TAVI. (See label and explanation of the symbols).

Contraindications

- Application in the central nervous system
- Guidewires with outer diameter $>0.014"$ (0.35mm) should not be used in the coronary vasculature.

Since guidewires are used to enable the performance of the above-mentioned procedures, the contraindication for these procedures also apply for the guidewires.

Application

General

The application of the guidewires is exclusively restricted to medical specialists with experience in performing the above-mentioned procedures.

The guidewires are intended for patients who are qualified to undergo the above-mentioned procedures. The decision whether a patient is qualified to undergo the procedures shall be made by the physician. The guidewires are intended for single use. The duration of use must not exceed 24 hours

The interventional combination device has to be prepared according to manufacturer's instructions of the device.

Summary of Safety and Clinical Performance

The Summary of Safety and Clinical performance can be found under the following link in EUDAMED: <https://ec.europa.eu/tools/eudamed>. The following Basic UDI-DI can be used to search for the SSCP in EUDAMED: 42556056gwclIIEH.

Warning Notices

- If used improperly, guide wires can lead to perforation of tissue. The consequences of perforation can be severe and can lead to death of the patient. Tissue perforation primarily occurs due to kinking of the guide wire. When applying the guide wires an insertion technique must be used that avoids kinking of the guide wire. An insertion applicator has to be used. Furthermore, there are known cases of overinsertion of guide wires; this is generally considered as a serious risk. Overinsertion of guide wires can lead to injury and/or perforation of the treated vessel.
- **ATTENTION:** Guide wires contain ferromagnetic materials and must not be used under MR conditions as they may move and/or heat up in magnetic fields. This can lead to catastrophic complications. Only guidewires labelled according can be used in the MRI
- The guide wires marked as „MR conditional“ (see symbol description) were tested for safety in applications under MR conditions and were

classified as MR conditional. They can be used under conditions specified below under MRT conditions:

- o Static magnetic fields of 1.5 T and 3 T
- o Maximum spatial field gradient of 12,900 G/cm (129 T/m)
- o Maximum force product of 229 T²/m
- o Theoretically expected maximum whole-body average (WBA) specific absorption rate (SAR) of 2 W/kg (normal operation mode)

- Under the conditions described above, after 15 minutes of continuous scanning, a maximum temperature rise is to be expected of less than
 - o 1,4 °C (2 W/kg, 1.5 Tesla) RF-dependent temperature increase with a background temperature rise of ≈ 1.3 °C (2 W/kg, 1.5 Tesla)
 - o 2,0 °C (2 W/kg, 3 Tesla) RF-dependent temperature increase with a background temperature rise of $\approx 2,1$ °C (2 W/kg, 3 Tesla)
- In non-clinical tests, the image artifacts caused by the guide wire expanded by approximately 58,4 mm with a gradient echo pulse sequence and a 3 Tesla system.
- Repeated advancing and retracting of the guidewire should be avoided. This may lead to the formation of knots and entrapment of the guidewire. The user should be familiar with appropriate techniques in order to be able to control serious complications in such cases.
- If coated guide wires are used with other medical devices (e.g. catheter), special attention has to be paid to an appropriate combination to avoid damages of the coating. The important dimensions of the guide wire are indicated on the label.
- Certain cases of entrapped or broken guide wires have been reported as well as possible techniques to retrieve the entrapped or broken fragments. The user should be familiar with appropriate techniques in order to be able to control serious consequences in such cases.
- Like almost every product that is introduced into blood vessels, guide wires may lead to thrombosis or serious bleeding at the puncture point. The user is responsible for taking adequate measures to prevent these complications during the application of the guide wire.
- The guide wire should be manipulated slowly and carefully. The position and advancement of the guidewire tip should be controlled using imaging techniques such as fluoroscopy. Manipulation of the guidewire without controlling its position can lead to incorrect positioning of the guide-wire and injury to the patient.
- The guidewire must be held in position when changing or withdrawing the catheter in order to prevent injury to the patient.
- The markings on the guide wires are for orientation purposes only. They must not be used as a measuring aid.
- If any resistance is felt during use, the manipulation of the guide wire and/or the combination device used, such as catheter or endoscope, must be stopped immediately and the cause for the resistance must be determined by use of various visual methods. Otherwise, overexpansion, kinking, breaks or partial detachment of the outer polymer jacket or the coating material of the guide wire can occur. Remove, if re-quired, the guide wire and the combination devices used as a complete unit to avoid complications.
- The guide wire must not be withdrawn through a metallic cannulas or devices with sharp edges. The sharp edges of these devices can damage or shear the coating or polymer jacket.
- Extreme caution is required when using a device that emits energy (laser, pressure, ultrasound, etc.). The guide wire has to be retracted into a position where it cannot be damaged. Direct contact may cause damage to the wire and/or sever the wire.
- If the guide wire is inserted into the right atrium, there is a risk of cardiac arrhythmia or injury. It is recommended to monitor the patient during the procedure with ECG.
- If the guidewire is inserted into the left ventricle, there is a risk that it will penetrate the left atrium or tangle with the papillary muscle of the mitral valve. Check the anatomy of the left ventricle and select an appropriate guidewire.

• Guidewires made of stainless steel may contain cobalt which is classified as a Carcinogenic category 1B, Mutagenic category 2 toxic to re-production category 1B. However, based on the current data, there is no evidence that the cobalt in stainless steel results in a toxic, carcinogenic or reproductive toxic effect.

Preventive Measures

- Guidewires are sensitive devices which should be used very carefully. The device should only be used by experienced medical specialists. Furthermore, the application of the device requires a detailed understanding of the technical principles as well as the clinical applications and risks associated with the use of guide wires in order to avoid any damage to the guide wire and to prevent patient hazards.
- Before and if possible during its use the guide wire should be checked visually for damages, misalignment, kinks or abrasion of the outer polymer jacket or the coating and for other deformations. Damaged, misaligned, kinked and otherwise deformed guide wires must not be used any longer and must be discarded due to potential risk of tissue perforation. Likewise, guide wires with abrasion of the outer polymer jacket or the coating must not be used since this may cause detachment of the fragments.
- Prior to use check the compatibility of the guidewire with the possible combination devices
- If metal insertion aids are used during application, we recommend removing these metal insertion aids immediately after inserting the guide wire and replacing them with a plastic sheath, plastic canula or similar insertion aids.
- Due to variations in stiffness, inner diameter or shape of certain catheter tips, abrasions of the outer polymer jacket or the coating may occur when manipulating. If any resistance is felt during introduction of the catheter, such catheters should not be used.
- Confirmation of the guidewire position by ultrasound may reduce the risk of vascular perforation.
- Guidewires with complete or partial hydrophilic coating must not come into contact with alcohol, antiseptic substances or substances containing solvents.
- Guidewires with soft distal tip are recommended for anatomies with extreme dilations or sever angulation to avoid injuries.
- Cases of guidewire migrating into the patient have been reported in the literature. To prevent this from happening hold on to the guidewire at all time during the entire procedure
- Cases of an entire guidewire or broken parts of guidewires remaining in the vasculature system have been reported in the literature. Always inspect the guidewire for complete removal at the end of the procedure.

Complications

Complications that may result from the use of the guidewire include, but not limited to:

- Allergic reaction
- Bleeding at the access site
- Infection
- Vessel trauma (Dissection, Perforation, injury)
- Fracture/breakage of the guidewire
- Entrapment of guidewire
- Damage to coating or the polymer jacket

Sterility

The guidewires are supplied sterile. The sterilization method is EO gas.

Handling

The guidewire is delivered in a sterile state and must be used in an aseptic environment. The guidewire should immediately be presented to the aseptic operation area after opening the packaging

The device is only sterile when the packaging is not damaged and not opened as well as when the expiry date has not been exceeded. A guidewire whose packaging is damaged or its expiry date is exceeded must not be used.

The guidewires are intended for single use and should not be reprocessed for re-use. Reprocessing can compromise the sterility and/or the functionality of the device

Storage

The guidewire must be stored in a dry place away from direct light.

Disposal


After the procedure the guidewire must be disposed of in a safely and orderly manner according to the applicable legal requirements for the disposal of medical devices. The MR guide wires must not be shredded.

Report of serious incidences

Any serious incidence that occurs in relation to the devices should be reported to the manufacturer and to the responsible competent authority.

Annotation

Symbols used on the packaging and labelling of these devices have the following meaning:

Symbol	Denotation	Symbol	Denotation
	Date of manufacture		Contains hazardous substances
	Manufacturer		Distributor
	Use by date		Single sterile barrier system
	Batch code		Medical device
	Catalogue number		Refer to instruction manual/booklet
	Sterilized using ethylene oxide		MR unsafe. Do not use this device under MR conditions
	Keep dry		MR Conditional *
	Keep away from sunlight		Non-Pyrogenic
	Do not use if package is damaged and consult instruction for use		TAVI application
	Do not reuse		Authorized representative
	Do not re-sterilize		

* Can be used in combination with 1.5 Tesla and 3 Tesla devices under the conditions described above.

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EPflex Feinwerktechnik GmbH