EPflex

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ENGLISH

Note Carefully read and understand the entire Instruction for use prior to using the guidewire. Failure to do so could result in complication.

Intended purpose

The guidewires are intended to facilitate the introduction and placement of diagnostic or therapeutic devices in peripheral vasculature or hollow organs of the human body during endoscopic or interventional procedures.

Indications

The guidewires are indicated for use during the following clinical procedures:

- Percutaneous Transluminal Angioplasty
- Ureterorenoscopy
- Percutaneous Nephrolithotomy
- Percutaneous Nephrostomy
- Placement of devices in the urinary tract e.g. ureter catheter
- Endoscopic Retrograde Cholangio-Pancreatography
- Placement of drainage catheters in different hollow organs of the body involving the use of Seldiger technic

ATTENTION: A guidewire may only be used for the procedure for which it is labelled. See label and symbol description.

Contraindications

- Application involving the direct contact central nervous system
- Application involving the direct contact with the heart or the central circulatory system
- Guidewires with hydrophilic coating along the entire length may not be used in the peripheral vasculature.

In general, since guidewires for Peripheral Application are used to enable the performance of the above-mentioned procedures, the contraindication for these procedures also apply for the guidewires.

Application

General

The application of the guidewires is exclusively restricted to medical specialists with experience in performing the above-mentioned procedures.

The guidewires are intended for patients who are qualified to undergo the above-mentioned procedures. The decision whether a patient is qualified to undergo the procedures shall be made by the physician. The guidewires are intended for single use. The duration of use must not exceed 24 hours

The interventional combination device has to be prepared according to manufacturer's instructions of the device.

Guide Wires with Hydrophilic Coating

Before removing the fully hydrophilic coated guidewire from the dispenser, physiological saline solution has to be injected into the hub end of the dispenser in order to wet the complete hydrophilic surface of the guide wire.

After the injection of the saline solution, the guidewire has to be withdrawn carefully from the dispenser. If the guidewire cannot be withdrawn easily from the dispenser more saline solution has to be injected into the dispenser.

After removal the guidewire must not be reinserted into the dispenser. The wetted guidewire must be used quickly and should not be allowed to dry, as this may result in it sticking to the catheters placed over it.

Warning Notices

- If used improperly, guidewires can lead to perforation of tissue. The consequences of perforation can be severe and can lead to death of the patient. Tissue perforation primarily occurs due to kinking of the guide wire. When applying the guidewires an insertion technique

must be used that avoids kinking of the guidewire. An insertion applicator has to be used. Furthermore, there are known cases of overinsertion of guidewires; this is generally considered as a serious risk. Overinsertion of guidewires can lead to injury and/or perforation of the treated vessel.

- **ATTENTION: These guidewires must not be used in MRI!** Guidewires contain ferromagnetic materials and must not be used under MR conditions as they may move and/or heat up in magnetic fields. This can lead to catastrophic complications.
- Repeated advancing and retracting of the guidewire should be avoided. This may lead to the formation of knots and entrapment of the guidewire. The user should be familiar with appropriate techniques in order to be able to control serious complications in such cases.
- If coated guidewires are used with other medical devices (e.g. catheter), special attention has to be paid to an appropriate combination to avoid damages of the coating. The important dimensions of the guidewire are indicated on the label.
- Certain cases of entrapped or broken guidewires have been reported as well as possible techniques to retrieve the entrapped or broken fragments. The user should be familiar with appropriate techniques in order to be able to control serious consequences in such cases.
- Like almost every product that is introduced into blood vessels, guidewires may lead to thrombosis or serious bleeding at the puncture point. The user is responsible for taking adequate measures to prevent these complications during the application of the guidewire.
- The guidewire should be manipulated slowly and carefully. The position and advancement of the guidewire tip should be controlled using imaging techniques such as fluoroscopy. Manipulation of the guidewire without controlling its position can lead to incorrect positioning of the guidewire and injury to the patient.
- The guidewire must be held in position when changing or withdrawing the catheter in order to prevent injury to the patient.
- The markings on the guidewires are for orientation purposes only. They must not be used as a measuring aid.
- If any resistance is felt during use, the manipulation of the guidewire and/or the combination device used, such as catheter or endoscope, must be stopped immediately and the cause for the resistance must be determined by use of various visual methods. Otherwise, overexpansion, kinking, breaks or partial detachment of the outer polymer jacket or the coating material of the guidewire can occur. Remove, if required, the guidewire and the combination devices used as a complete unit to avoid complications.
- The guidewire must not be withdrawn through a metallic cannulas or devices with sharp edges. The sharp edges of these devices can damage or shear the coating or polymer jacket.
- Extreme caution is required when using a device that emits energy (laser, pressure, ultrasound, etc.). The guidewire has to be retracted into a position where it cannot be damaged. Direct contact may cause damage to the guidewire and/or sever the guidewire.
- Guidewires made of stainless steel may contain cobalt which is classified as a Carcinogenic category 1B, Mutagenic category 2 toxic to reproduction category 1B. However, based on the current data, there is no evidence that the cobalt in stainless steel results in a toxic, carcinogenic or reproductive toxic effect.

Preventive Measures

- Guide wires are sensitive devices which should be used very carefully. The device should only be used by experienced medical specialists. Furthermore, the application of the device requires a detailed understanding of the technical principles as well as the clinical applications and risks associated with the use of guide wires in order to avoid any damage to the guide wire and to prevent patient hazards.
- Before and if possible during its use the guide wire should be checked visually for damages, misalignment, kinks or abrasion of the outer polymer jacket or the coating and for other deformations.

Damaged, misaligned, kinked and otherwise deformed guide wires must not be used any longer and must be discarded due to potential risk of tissue perforation. Likewise, guide wires with abrasion of the outer polymer jacket or the coating must not be used since this may cause detachment of the fragments.

- Prior to use check the compatibility of the guidewire with the possible combination devices
- If metal insertion aids are used during application, we recommend removing these metal insertion aids immediately after inserting the guide wire and replacing them with a plastic sheath, plastic canula or similar insertion aids.
- Due to variations in stiffness, inner diameter or shape of certain catheter tips, abrasions of the outer polymer jacket or the coating may occur during manipulation. If any resistance is felt during introduction of the catheter, such catheters should not be used.
- Guide wires with complete or partial hydrophilic coating must not come into contact with alcohol, antiseptic substances or substances containing solvents.
- Cases of guidewire migrating into the patient have been reported in the literature. To prevent this from happening hold on to the guidewire at all time during the entire procedure
- Cases of an entire guidewire or broken parts of guidewires remaining in the vasculature system have been reported in the literature. Always inspect the guidewire for complete removal at the end of the procedure.

Complications

Complications that may result from the use of the guidewire include, but not limited to:

- Allergic reaction
- Bleeding at the access site
- Infection
- Vessel trauma (Dissection, Perforation, injury)
- Fracture/breakage of the guidewire
- Entrapment of guidewire
- Damage to coating or the polymer jacket

Sterility

The guidewires are supplied sterile. The sterilization method is EO gas.

Handling

The guidewire is delivered in a sterile state and must be used in an aseptic environment. The guidewire should immediately be presented to the aseptic operation area after opening the packaging.

The device is only sterile when the packaging is not damaged and not opened as well as when the expiry date has not been exceeded. A guidewire whose packaging is damaged or its expiry date is exceeded must not be used.

The guidewires are intended for single use and should not be reprocessed for re-use. Reprocessing can compromise the sterility and/or the functionality of the device

Storage

The guidewire must be stored in a dry place away from direct light.

Disposal

After the procedure the guidewire must be disposed of in a safely and orderly manner according to the applicable legal requirements for the disposal of medical devices.

Report of serious incidences

Any serious incidence that occurs in relation to the devices should be reported to the manufacturer and to the responsible competent authority.

Annotation

Symbols used on the packaging and labelling of these devices have the following meaning:

Symbol	Denotation	Symbol	Denotation
	Date of manufacture		Distributor
	Manufacturer		Single sterile barrier system
	Use by date		Medical device
	Batch code		Refer to instruction manual/booklet
	Catalogue number		MR unsafe Do not use this device under MR conditions
	Sterilized using ethylene oxide		Non-Pyrogenic
	Keep dry		Urological application
	Keep away from sunlight		Gastrological application
	Do not use if package is damaged and consult instruction for use		Peripheral vascular application
	Do not reuse		Authorized representative
	Do not re-sterilize		Contains hazardous substances

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EPflex Feinwerktechnik GmbH