EPflex

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## ENGLISH

**Note:** *Carefully read and understand the entire instruction for use prior to using the stone retrieval device. Failure to do so could result in complications*

### Device Description

The stone retrieval devices consist of four main parts: an expandable stone dislodger such as a basket or a grasper, a pull wire, a flexible sheath and a handle. The self-expanding dislodger at the distal end of the device is connected with the handle at the proximal end of the device with the pull wire housed in the flexible sheath. The handle is used to close and open the dislodger during use.

The stone retrieval devices are intended to be used through the working channels of an endoscope.

### Intended Purpose

Stone Retrieval Devices are intended for the endoscopic removal of stones from the urogenital tract or the common bile duct during retrograde interventions.

### Indications

The SRD are indicated for medical conditions requiring the removal of stones using the following procedures:

- Ureterorenoscopy (URS)
- Endoscopic Retrograde Cholangio-Pancreatography (ERCP)

**ATTENTION:** A stone Retrieval Device may only be used for the procedures for which it is labeled. See label and symbol description.

### Contraindications

The SRD must not be used for:

- Percutaneous Nephrolithotomy (PCNL)
- Mechanical stone crushing (lithotripsy)

The contraindication for the above-mentioned endoscopic procedures apply.

### Application

#### Intended User

The application of the devices is exclusively restricted to medical specialists with experience in performing the above-mentioned procedures.

#### Duration of use

The stone retrieval devices are intended for single use. The duration of use must not exceed one hour.

### Application Instruction

#### Prior to Use

1. Make sure the sterile packaging is not damaged. Do not use the device if the sterile packaging is damaged.
2. Present the device to a sterile field.
3. Inspect the device for any visible damages.
4. Test the device by opening and closing the dislodger several times before introducing it into the patient.

#### Insertion

5. Ensure that the dislodger is closed before inserting the stone retrieval device into the working channel of an endoscope.
6. Insert the stone retrieval device with the closed dislodger into the working channel of the endoscope.
7. Carefully advance the stone retrieval device under visual control through the working channel of the endoscope until it can be seen exiting the distal end of the endoscope. During the advancement ensure that the dislodger remains closed.

#### Capture

8. Advance the stone retrieval device with the closed dislodger under visual control towards the stone. In the case of baskets, the dislodger is positioned above the stone, in the case of graspers the

dislodger is positioned in front of the stone.

9. Once the dislodger is in the right position move the slider forward to open the dislodger.
10. Maneuver the stone retrieval device under visual control to place the stone in the opened dislodger.
11. Once the stone is in the opened dislodger close the dislodger slowly under visual control by moving the slider backwards until the stone is gripped tightly.

#### Removal

12. Once the stone is captured carefully withdraw the stone retrieval device under visual control, thereby removing the stone from the body.

### Warnings

- If used improperly, stone retrieval devices can cause the perforation of tissue, in particular if the stones are lodged on the vessel wall. The stone bed is then frequently very fragile.
- Excessive mechanical forces can tear or kink the stone retrieval device.
- Kinks on the shaft impede the mechanical functionality of the stone retrieval device.
- If the stone retrieval device is severely bent, the deployment or retraction of the dislodger becomes difficult or impossible. In this case, the en-tire stone retrieval device should be pushed forwards or backwards a few millimeters.
- The use of contrast media can lead to adhesions that can limit the functionality of the stone retrieval device.
- When using products emitting energy (laser, pressure, ultrasound....) extreme caution should be exercised and the instrument retracted to a position in which it cannot be damaged. Direct contact with these devices can damage and/or destroy the stone retrieval device.
- Some stones may be too large to be removed with the stone retrieval device through the endoscope because the stone could get stuck in the working channel of the endoscope during removal.
- The stone retrieval devices must not be used under MR-field.
- Stone retrieval devices may not be used for mechanical stone crushing (lithotripsy).
- Do not treat products with disinfectant; no data concerning material changes is available to EPflox.
- Stone Retrieval Devices with stainless steel components may contain cobalt, which is classified as carcinogen category 1B, mutagenic category 2 and toxic to reproduction category 1B. However, based on current data, there is no evidence that the cobalt in stainless steel results in a toxic, carcinogenic or reproductive toxic effect.

### Precautionary Measures

- Stone retrieval devices are sensitive instruments and should be handled with extreme care. Use of the product presupposes a detailed understanding of the technical principles, clinical application options and risks associated with using stone retrieval devices, so as to prevent damage to the instrument and avoid risk to the patient. The product may therefore only be used by medical professionals who are experienced in the handling of this medical product and the relevant surgical techniques.
- When withdrawing the stone retrieval device, the shaft must be kept as straight as possible. Withdrawing the stone retrieval device at excessive angle or curvature can damage the shaft.
- The size of the stone should be determined using imaging techniques before the intervention so as to prevent entrapment. Do not use the stone retrieval device if the stone is too large.
- The size of the stone should be determined using imaging techniques before the intervention to ensure a right device is chosen. Do not use the stone retrieval device if the stone is too large.
- As soon as resistance is felt during extraction with the stone retrieval device, the stone retrieval device must be loosened and the stone disintegrated. If this is not successful, do not on any account use force for the extraction; there is a high risk of the vessel wall tearing off. In this case, the stone in the instrument must be subjected to

lithotripsy. After the lithotripsy the stone retrieval device can be removed.

- Before use, the function (opening/closing) of the instrument should always be checked. If possible, the stone retrieval device should be inspected carefully for damage, distortion, kinks and other deformations during the procedure. Damaged, distorted, kinked, and otherwise deformed stone retrieval devices may not be used and must be disposed of.
- The compatibility of the stone retrieval device with the working channel of the endoscope should be checked before the intervention.
- Advance or remove the instrument carefully in order to prevent injury to the mucous membrane or the urinary tract / bile duct system (trauma, perforation). The advancement or removal should always be done under visual control to avoid injuries to the patient.
- Extreme caution should be exercised in case of simultaneous laser lithotripsy application, as the nitinol instruments can be damaged very easily by the laser. Do not on any account point the laser directly into the instrument.

### Possible complications


The following complications are possible when using stone retrieval devices for stone removal:

- Entrapment of large stones
- Inability to disengage the dislodger from irretrievable stones requiring the application of other interventions
- Tissue perforation
- Breakage of the stone retrieval device
- Infection
- Non-retrievable stones.


### Disassembling of the EPflex System handle. Only for Stone Retrieval Devices with EPsys handle.

In case the dislodger needs to be separated from the handle during the procedure, perform the following steps


Push the release button on the underside of the handle

Lift the sliding unit from the handle

To disconnect the pull wire from the sliding unit, insert the screwing-element integrated in the Handle in the eccentric screw and turn it anti-clockwise

Pull the thumb-slide out of the unit

Remove the sheath by pulling the part connected to it backwards. The dislodger and the pull wire remain in the body

### Sterility

The stone retrieval devices are supplied sterile. The sterilization method is EO gas.

### Handling

The stone retrieval devices must be used in an aseptic environment. The stone retrieval devices should immediately be presented to the aseptic operation area after opening the packaging.

The sterile of the stone retrieval devices only guaranteed when the sterile packaging is not damaged and not opened as well as when the expiry date has not been exceeded. A stone retrieval device whose sterile packaging is damaged or its expiry date is exceeded must not be used.

The stone retrieval devices is intended for single use and should not be reprocessed for re-use. Reprocessing can compromise the sterility and/or the functionality of the device.

### Storage

The stone retrieval devices must be stored in a dry place away from direct light.

### Disposal

After the procedure the stone retrieval devices and the packaging materials must be disposed of in a safely and orderly manner according to the applicable hospital and/or legal requirements for the disposal of medical devices.

### Report of serious incidence

Any serious incidence that occurs in relation to the devices should be reported to the manufacturer and to the responsible competent authority.

### Annotation

Symbols used on the packaging and labelling of these devices have the following meaning:

Symbol	Denotation	Symbol	Denotation
	Date of manufacture		Distributor
	Manufacturer		Single sterile barrier system
	Use by date		Medical device
	Batch code		Refer to instruction manual/ booklet
	Catalogue number		MR unsafe Do not use this device under MR conditions
	Sterilized using ethylene oxide		Non-Pyrogenic
	Keep dry		Urological application
	Keep away from sunlight		Gastrological application
	Do not use if package is damaged and consult instruction for use		Authorized representative
	Do not reuse		Contains hazardous substances
	Do not re-sterilize		

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EPflex Feinwerktechnik GmbH